

**Product Name :**  
Nauga Lab Slide Warming Table

**Product Code :**  
NLAB-LABORATORYEM42003



**Description :**

Nauga Lab Slide Warming Table

**Technical Specification :**

Slide Warming Table

A versatile equipment to heat slides with just microtome sections to attain maximum flatness before Staining.

Temperature range ambient temperature is thermostatically controlled.

The body is made out of thick mild steel duly finished in stoving enamel and top is made of highly polished copper or stainless steel sheet.

Complete with thermostatic control knob, indicator lamp, plug & cord etc.

To work on 220 V AC 50 Hz single phase.

Size In cm Size in Inches

30x25 (12"x10")

45x30 (18"x12")

60x15cm (24"x6")

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003