

Product Name : Naugra Lab Standard X-Block	Product Code : NLAB-LABORATORYEM77017
Description : Naugra Lab Standard X-Block	
Technical Specification : Standard X-Block : Machined from solid Brass rod, richly chromium plated to hold 9 mm & 6 mm rods. Standard Large X-Block: Machined from solid Brass rod, richly chromium plated to hold 12.5mm & 9mm rods	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003