

<b>Product Name :</b> Naugra Lab Standard X-Block	<b>Product Code :</b> NLAB-LABORATORYEM77017
--	---

--

**Description :**

Naugra Lab Standard X-Block

**Technical Specification :**

Standard X-Block :

Machined from solid Brass rod, richly chromium plated to hold 9 mm & 6 mm rods.

Standard Large X-Block:

Machined from solid Brass rod, richly chromium plated to hold 12.5mm & 9mm rods

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003