

Product Name :
Nauga Lab Student Induction Coil

Product Code :
NLAB-LABORATORYEM21022



Description :

Nauga Lab Student Induction Coil

Technical Specification :

Student Induction Coil

On teakwood base with graduated scale.

Primary coil with 300 turns and secondary with 6500 turns.

Interrupter working directly off the iron core.

Complete with heavy duty contacts.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003