


Product Name :
Naugra Lab Student Kymographs Digital

Product Code :
NLAB-LABORATORYEM77003



Description :

Naugra Lab Student Kymographs Digital

Technical Specification :

Student
Kymographs (Recording Drum-Digital)

Digital Microprocessor controlled
1st time in the world of Pharmacology Digital kymograph/Sherrington Drum
Designed and Manufactured by Microteknik.

Detail Description

Digital kymographs:- 16 x 2 LCD
Digital Display with Back Light simultaneous display of all functions on the
display and selection through display.

Digital Speed Selection

Highly Accurate Speed : E/C
System Speeds Available -0.12, 0.25, 0.50, 0.75, 1.00, 1.25, mm/Sec With
Accuracy of +/-0.01%

Very Low Power Consumption (12
volts 1 Amp.)

Operated with Stepper Motor which
is highly accurate and run by Microprocessor based program.

Special Chrome Plated and Powder
Coated Body:Extremely Resistant to Corrosion. Self Lubricating Gear.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003