

Product Name :
Naugra Lab Table Top Autoclave

Product Code :
NLAB-LABORATORYEM47007



Description :

Naugra Lab Table Top Autoclave

Technical Specification :

Table Top Autoclave

Chamber capacity: 20L Automatic control and water inlet Auto-dry mode

Overall dimensions: 624x430x621mm (WxHxD)

Power Consumption: 220VAC, 50/60Hz, 2000W Sterilization

Temperature Pressure: 121oC(249.8oF)/1.1kg/cm2/15.64(lb/in2) 134oC(273.2oF)/2.1kg/cm2/29.85(lb/in2)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003