

**Product Name :**  
Nauga Lab Tablet Hardness Tester

**Product Code :**  
NLAB-LABORATORYEM24003



**Description :**

Nauga Lab Tablet Hardness Tester

**Technical Specification :**

Tablet Hardness Tester (Monsanto type) : Useful to take quick readings from 0-20 Kg pressure for hardness. The equipment is made of brass and richly chrome plated. The index and division of reading is 0.1 Kg. Completed in box

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003