

Product Name :
Nauga Lab Tablet Hardness Tester

Product Code :
NLAB-LABORATORYEM24003



Description :

Nauga Lab Tablet Hardness Tester

Technical Specification :

Tablet Hardness Tester (Monsanto type) : Useful to take quick readings from 0-20 Kg pressure for hardness. The equipment is made of brass and richly chrome plated. The index and division of reading is 0.1 Kg. Completed in box

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003