

Product Name :
Naugra Lab Tablet Making Machine Motorised

Product Code :
NLAB-LABORATORYEM23003



Description :

Naugra Lab Tablet Making Machine Motorised

Technical Specification :

Hammer type and provided with suitable motor, workable on 220V AC 50Hz single phase.

Optional : Spare punch & die of chrome steel of 6, 9 & 12 mm.

Note : Ask for rotary tablet machine

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003