

Product Name :
Naugra Lab Tambours

Product Code :
NLAB-LABORATORYEM77011

Description :

Naugra Lab Tambours

Technical Specification :

Tambours

Differnt types of Tambours are available with following features:-

Student Tambour:- Size-40 mm tie on diaphragm, Stainless Steel capillary lever.

Large Tambour:- 75mm diameter diaphragm for large displacement of air, Stainless Steel capillary lever.

Marey's Tambour:- Diameter 22 mm or 28 mm with Stainless Steel capillary lever

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003