

Product Name :
Naugra Lab Temperature and Precision Control Water Bath

Product Code :
NLAB-LABORATORYEM46021



Description :

Naugra Lab Temperature and Precision Control Water Bath

Technical Specification :

Temperature and Precision Control Water Bath

Specifications of Water bath but temperature controlled by solid state controller with stirrer with S.S. rod & blade to make temp. of bath all-round uniform and constant.

Chamber of stainless steel sheet.

Temperature Range:- 50° C above ambient to 100° C.

Chamber Size:- 430 x 240 x 265 mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003