

**Product Name :**  
Naugra Lab Tissue Culture Rack

**Product Code :**  
NLAB-LABORATORYEM6604



**Description :**

Naugra Lab Tissue Culture Rack

**Technical Specification :**

**Tissue Culture Rack**

Tissue Culture Rack made out of mild steel pipe frame.

Distance between shelves to shelves is 16".

Size of each shelf is 48" x 18".

Platform surface will be of 3mm thick hylem Sheet.

Four 40-watt fluorescent light (With Electronic Ballast).

Minimum heating effect/heat dissipation.

No Flickering effect.

Made out of mild steel pipe frame.

Distance between shelves to shelves is 16".

Size of each shelf is 48" x 18".

Platform surface will be of 3mm thick hylem Sheet.

Four 40-watt fluorescent light (With Electronic Ballast).

Minimum heating effect/heat dissipation.

No Flickering effect.

Energy saver with individual switching arrangement in each shelf

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003