

**Product Name :**  
Naugra Lab Tissue Floatation Bath

**Product Code :**  
NLAB-LABORATORYEM46014



**Description :**

Naugra Lab Tissue Floatation Bath

**Technical Specification :**

Tissue Floatation Bath

Tissue floatation water bath is most ideal for Histological, chemical, clinical and bacteriological laboratories used to stretch and dry cut tissues.

Special bath reliably fulfills everyday task in laboratories.

Its exact temperature regulation ensures evenly stretched cuttings which are clearly visible in the anodised bath.

The cuttings are carefully dried on the warmed rim of the bath.

Offers highest level of performance, flexibility and control for most demanding applications.

Exact temperature constancy of  $\pm 0.5^{\circ}\text{C}$ .

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003