

Product Name :
Nauga Lab Tissue Homogeniser

Product Code :
NLAB-LABORATORYEM50001



Description :

Nauga Lab Tissue Homogeniser

Technical Specification :

Tissue Homogeniser

Tested on varied qualitative measures, this homogenizer is manufactured by utilizing quality approved components under the surveillance of our adroit experts.

The offered homogenizer is mainly used in research labs for the homogenization of tissue, plant, food, soil and various other materials.

Further, the provided Laboratory Tissue Homogenizer can be availed from us at highly competitive rates.

Features:-

Fitted with a holder

Sturdy design

Optimal strength

Suitable for tissue
grinding

Homogenizing Vessel
25 ml
50 ml
100ml

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003