

Product Name :
Nauga Lab Tray Drying Oven

Product Code :
NLAB-LABORATORYEM40003



Description :

Nauga Lab Tray Drying Oven

Technical Specification :

Tray Drying Oven

Suitable for various industrial, laboratory and R&D applications.

Heavy Forced Air Circulation by Strong and Reliable axial flow fans for cross flow circulation.

Voltage and Ampere meter equipped for monitoring safe operating of the machine.

Digital PID Control System enables superior temperature accuracy.

Digital Timer function with Buzzer Alarm (Option).

Over Temperature & Over Current Protection.

Adjustable Damper to exhaust moist air or fumes/vapours which in turn enhances thermal efficiency.

Heating is through special grade tubular air heaters placed in the circulating air path.

Sturdy triple walled unit built with suitable mild steel sections and heavy gauge mild steel sheet

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003