


Product Name :
Naugra Lab True Double Beam Atomic Absorption

Product Code :
NLAB-LABORATORYEM76007



Description :

Naugra Lab True Double Beam Atomic Absorption

Technical Specification :

True Double Beam Atomic Absorption Spectrophotometer

Salient Features & Benefits

High Resolution Touch Screen LCD Color Graphical Display Supports USB Pen Drive To Store / Retrieve MASS Spectral Data.

BUILT-IN PC Control with Advanced capability for Data Processing, Processing, Storage, Retrieval, Interpretation, etc.

Background correction with D2 lamp source.

Auto Burner Height Optimization, and maximum ABS adjustment of burner through membrane keyboard .

Auto Flame control.

Integral Microcontroller based auto gas box with all safety features.

Supports USB Printing with PCL Class Printers.

Applications For analysis of various metallic constituents in:-

Soil

Fertilizer

Micro Nutrients

Plants

Pesticide

Foods

Water

Biological

Matter

Metallurgy

Plating Coal

Petrochemicals

Glass

Cement

Steel

Others

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003