

**Product Name :**  
Nauga Lab Universal Moisture Meter Digital

**Product Code :**  
NLAB-LABORATORYEM51003



**Description :**

Nauga Lab Universal Moisture Meter Digital

**Technical Specification :**

**Universal Moisture Meter Digital:-**

The most versatile instrument for the measurement of Moisture content in Grains, Seeds and Pulses.

Operates on 220 V 50Hz and Battery.

The instrument is made to operate very easily and takes less than a minute to complete the test.

Rugged in construction, aluminum casted and the electrical unit is contained in.

The bevel gear mechanism  
is of best steel.

The whole instrument is  
protected in wooden box.

Note:- The instrument is  
approved by food corporation of India Central Ware Housing Corporation, Food  
Supply Department & Other Agriculture Institutions.

#### Specification:-

Measuring Range	0-100%
Accuracy	±0.2%
Dimension	32 x 17 cm
Weight	1050 Kg. without wooden cover.

Power Supply Cells or 220 V 50 Hz
---

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003