


Product Name :
Nauga Lab Universal Moisture Meter Digital

Product Code :
NLAB-LABORATORYEM51003



Description :

Nauga Lab Universal Moisture Meter Digital

Technical Specification :

Universal Moisture Meter Digital:-

The most versatile instrument for the measurement of Moisture content in Grains, Seeds and Pulses.

Operates on 220 V 50Hz and Battery.

The instrument is made to operate very easily and takes less than a minute to complete the test.

Rugged in construction, aluminum casted and the electrical unit is contained in.

The bevel gear mechanism is of best steel.

The whole instrument is protected in wooden box.

Note:- The instrument is approved by food corporation of India Central Ware Housing Corporation, Food Supply Department & Other Agriculture Institutions.

Specification:-

Measuring Range 0-100%	
Accuracy ±0.2%	
Dimension 32 x 17 cm	
Weight without wooden cover.	

Power Battery Cells or 220 V 50 Hz
--

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003