

Product Name :
Naugra Lab Vacuum Pump

Product Code :
NLAB-LABORATORYEM53001



Description :

Naugra Lab Vacuum Pump

Technical Specification :

Vacuum Pump

Vacuum Pump Motor driven, oil sealed pump mounted on heavy base.

It is driven by V belt at speed up-to 400 RPM.

Choice of single and double stage available creating Ultimate vacuum of 0.01 mm of Hg column and 0.0 mm of Hg column respectively.

The pump is to work on 220V AC 50 Hz single phase.

Single Stage

Motor (HP)	
0.375	
0.50	
0.75	
1.50	

Double Stage	
Motor (HP)	
0.125	

0.125	
0.45	
0.45	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003