

Product Name :
Naugra Lab Variable Interrupter

Product Code :
NLAB-LABORATORYEM21023

Description :

Naugra Lab Variable Interrupter

Technical Specification :

Variable Interrupter

We also provide Tetanus set with Complete unit on two stands.

Mercury cup and magnet.

Other available recorders:-

Piston recorder

Out-flow (Float) recorder

Condons drop recorder

Inchleys drop recorder

Mechanical Strhomour

Jaquet chronometer

Jackson Enterograph

Oncometer kidney

Oncometer Heart

Hair Aesthesiometer.

Compass Aesthesiometer

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003