

Product Name :
Nauga Lab Vertical Laminar Air Flow Cabinet

Product Code :
NLAB-LABORATORYEM6601



Description :

Nauga Lab Vertical Laminar Air Flow Cabinet

Technical Specification :

Vertical Laminar Air Flow Cabinet

Specification:-

Vertical Laminar Air Flow Benches are used for various works that require maintenance of the cleanliness and purity like plant tissue culture.

Examination and experimentation of precision electronic parts, medical science, pharmacy, agriculture sector and the like.

Designed so as to meet the requirement of US Federal Standard 209 B providing particle free air to meet class 100 conditions.

Cabinet is fabricated of Wooden board laminated with sunmica OR Cold rolled steel duly powder coated OR Stainless Steel (as ordered).

Working table is made of Stainless Steel Sheet.

Side panels and front door are made out of thick transparent plexi glass.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003