

Product Name :
Naugra Lab Viscometric Water Bath

Product Code :
NLAB-LABORATORYEM46027



Description :

Naugra Lab Viscometric Water Bath

Technical Specification :

Viscometric Water Bath

Specification:-

To Check Kinematic Viscosity.

Double wall inner S.S., Outer G.I. Powder Coated.

Uniform Temperature with stirrer/circulatory pump.

Special grade immersion heater (Booster for initial boosting)

Glass viewing window.

Capacity (Viscosity holder)/6/9/12 or as per customer application.

Power 230V., A.C. 50Hz, mains.

Temperature ambient +50C to 900C, + 0.10C (0.010 on request)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003