

Product Name :
Naugra Lab Vortex Shaker

Product Code :
NLAB-LABORATORYEM39001



Description :

Naugra Lab Vortex Shaker

Technical Specification :

Vortex Shaker, Cyclo Mixer, Test Tube Shaker

Useful for accelerated mixing of solutions in test tubes, small flasks or centrifugal tubes etc.

The housing is 150x150 mm.

The speed is regulated by a speed control device.

Short eccentric stroke motor imparts vigorous agitation to the tubes or flasks.

Body is strong and finished in staving enamel paint.

The equipment is workable on 220V AC 50Hz single phase.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003