

Product Name :
Nauga Lab Water Bath Incubator Shaker Metabolic Shaker

Product Code :
NLAB-LABORATORYEM46013



Description :

Nauga Lab Water Bath Incubator Shaker Metabolic Shaker

Technical Specification :

Water Bath Incubator Shaker Metabolic Shaker

Water Bath Incubator Shaker has highest level of performance, flexibility and control for most demanding applications. Minimum Noise and Vibration is ensured by Specially-designed shaking mechanism.

Electronically Controlled Shaking Mechanism provides Quiet Reciprocating Motion and Precise Speed Control.

User's Self-compensation function to control difference between real temperature and indicated temperature.

Stainless Steel Bath for superior Durability & High Thermal efficiency.

Industrial grade energy efficient heaters reduces power bills drastically.

Double walled construction.

Space between the walls is filled with special grade mineral wool in order to avoid thermal losses.

Offers highest level of performance, flexibility and control for most demanding applications.

Minimum Noise and Vibration is ensured by Specially-designed shaking mechanism.

Electronically Controlled Shaking Mechanism provides Quiet Reciprocating Motion and Precise Speed Control.

User's Self-compensation function to control difference between real temperature and indicated temperature.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003