

**Product Name :**  
Naugra Lab Water Bath Paraffin

**Product Code :**  
NLAB-LABORATORYEM46006



**Description :**

Naugra Lab Water Bath Paraffin

**Technical Specification :**

Water Bath Paraffin

The Paraffin Section Flotation Bath is designed to assist with the handling of paraffin wax samples in histology and pathology laboratories.

It is essentially hot distilled water floating out bath that allows for the meticulous manipulation and location of sections onto glass slides.

The Paraffin Section Flotation bath has improved heat up times achieving a stable 45°C using 2.25 liters of water within 40 minutes.

The stability and accuracy have also been improved due to the change of control which now also allows the user to set the temperature as opposed to the arbitrary power, thereby, reducing the risk of overheating sections.

The Paraffin Section Flotation Bath has a PTFE black interior and powder-coated aluminium exterior, which makes it virtually maintenance-free, extremely easy to clean and allows for improved visibility of floating sections.

Technical Specification:-

Capacity
Temperature
Temperature Stability
Bath material base with PTFE coated interior
Dimensions (mm x mm x mm)
Weight kg
Heater Power

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003