

**Product Name :**  
Nauga Lab Water Bath Precision

**Product Code :**  
NLAB-LABORATORYEM46007



**Description :**

Nauga Lab Water Bath Precision

**Technical Specification :**

Water Bath Precision

Construction as per model but with glass windows on both opposite sides, for observation Temp.

Controlled by mercury contact thermometer with a sensitivity of 0.2°C.

Suitable for precision laboratory activities requiring constant temp.

With high degree of accuracy

Complete with stirrer with speed control regulator workable on 220V AC 50 Hz single phase.

Chamber approx.) Size (in mm)
425x250x225
450x375x375
{240x18"x18"} 450x450

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003