

Product Name :
Nauga Lab Water Bath Rectangular Double Walled

Product Code :
NLAB-LABORATORYEM46016



Description :

Nauga Lab Water Bath Rectangular Double Walled

Technical Specification :

Water Bath Rectangular Double Walled

Water Bath Rectangular Double Walled are Microprocessor Based PID Digital Temperature Indicator-cum-Controller maintains maximum temperature accuracy.

User's Self-compensation function to control difference between real temperature and indicated temperature.

Stainless Steel Bath for superior Durability & High Thermal efficiency.

Flat Stainless Steel Lid with Concentric rings prevents evaporation.

An industrial grade energy efficient heater reduces power bills drastically.

Double walled construction.

Microprocessor Based PID Digital Temperature Indicator-cum-Controller maintains maximum temperature accuracy.

User's Self-compensation function to control difference between real temperature and indicated temperature.

Stainless Steel Bath for superior Durability & High Thermal efficiency.

Flat Stainless Steel Lid with Concentric rings prevents evaporation.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003