

**Product Name :**  
Nauga Lab Water Bath Rectangular single wall

**Product Code :**  
NLAB-LABORATORYEM46005



**Description :**

Nauga Lab Water Bath Rectangular single wall

**Technical Specification :**

Water Bath Rectangular single wall

(Single wall) Body made out of thick stainless steel, top made out of chrome plated sheet concentric rings to cover 3"dia holes are provided.

Heating elements have an automatic cut-off device which operates when the bath gets dried

Equipment working on 220V AC 50Hz single phase

It is used to incubate samples in water at a constant temperature over a long period of time.

All water baths have a digital or an interface to allow users to set a desired temperature

Standard double wall construction

Inner being of stainless steel outer of mess duly staving paint finish and gap filled for temp.

Insulation with glass wool top cover with 75 mm (3") holes with concentric rings.

Size inside (Stainless Steel)
300x250x100mm
355x405x100mm

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003