

Product Name :
Naugra Lab Water Bath Rectangular

Product Code :
NLAB-LABORATORYEM46004



Description :

Naugra Lab Water Bath Rectangular

Technical Specification :

Water Bath Rectangular

It is used to incubate samples in water at a constant temperature over a long period of time.

All water baths have a digital or an interface to allow users to set a desired temperature.

Standard double wall construction.

Inner being of stainless steel outer of mess duly staving paint finish and gap filled for temp.

Insulation with glass wool top cover with 75 mm (3") holes with concentric rings.

Temp. Range 5oC above ambient to 95oC type element with accuracy of

Imported thermostat “EGO”/ “JUMO” German with the accuracy of $\pm 0.5^{\circ}\text{C}$.

Microprocessor PID Digital temp.

Indicator-Cum –Controller.

Digital Temp. Indicator.

Features

High operational fluency

Easy installation

High load bearing capacity

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003