

Product Name :
Nauga Lab Water Bath Seriolioica

Product Code :
NLAB-LABORATORYEM46028



Description :

Nauga Lab Water Bath Seriolioica

Technical Specification :

Water Bath Seriolioica

Emperature range from ambient temperature to 80degC.

Thermostatic control with accuracy of +/- 0.5degC.

Double walled , inside stainless steel and outside mild steel painted in epoxy powder coating .

Bath consists of two pilot lamps.

Temperature control knob and on/off switch.

To work on 220/230 volts A.C.

Supplied with or without stirring arrangement , without racks and thermometer.

Lid of water bath is made of stainless steel sheet

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003