

Product Name :
Nauga Lab Water Bath Serological

Product Code :
NLAB-LABORATORYEM46015



Description :

Nauga Lab Water Bath Serological

Technical Specification :

Water Bath Serological

Water Bath High Precision Ideal for Biotechnology, Clinical, Environmental, Medical, Petroleum, Pharmaceutical or Industrial Applications.

Offers highest level of performance, flexibility and control for most demanding applications.

Microprocessor Based PID Digital Temperature Indicator-cum-Controller maintains maximum temperature accuracy.

User's Self-compensation function to control difference between real temperature and indicated temperature.

Stainless Steel Bath for superior Durability & High Thermal efficiency.

Ideal for Biotechnology, Clinical, Environmental, Medical, Petroleum, Pharmaceutical or Industrial Applications.

Offers highest level of performance, flexibility and control for most demanding applications.

Water Bath High Precisions are microprocessor based PID Digital Temperature Indicator-cum-Controller maintains maximum temperature accuracy.

User's Self-compensation function to control difference between real temperature and indicated temperature.

Stainless Steel Bath for superior Durability & High Thermal efficiency.

Flat Stainless Steel Lid prevents evaporation and keeps constant temperature.

Industrial grade energy efficient heaters reduces power bills drastically.

Double walled construction.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003