

Product Name :
Nauga Lab Water Still Manasty Type

Product Code :
NLAB-LABORATORYEM31001



Description :

Nauga Lab Water Still Manasty Type

Technical Specification :

Water Still Manasty Type

Casted stainless Steel vessel is electrically operated for getting pyrogen-free distilled water automatically.

Automatically ejection device which operates automatically when water vessel gets dried.

Output Capacity (Approximately)

Output 3/4 Litres per hour

Output 5/6 Liters per hour

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003