

**Product Name :**  
Nauga Lab Water Still Table Pattern-Barnstead Type

**Product Code :**  
NLAB-LABORATORYEM71002



**Description :**

Nauga Lab Water Still Table Pattern-Barnstead Type

**Technical Specification :**

Water Still Table Pattern-Barnstead Type

Double walled construction with inner and outer chamber made of stainless steel.

The gap in between filled with glass wool.

Condenser & Top are also made of Stainless Steel.

SS-304 grade Argon Arc Welded.

Supplied with the rubber tubing coil of 5 meters.

Supplied with cord & plug to work on 220/230 volts AC supply.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003