

**Product Name :**  
Nauga Lab Water Still Wall Mounting

**Product Code :**  
NLAB-LABORATORYEM71003



**Description :**

Nauga Lab Water Still Wall Mounting

**Technical Specification :**

**Water Still Wall Mounting**

The unit is recommended for making pyrogen free distilled water for routine laboratory requirements.

Made out of Stainless Steel (SS-304 grade) sheet.

The cock fitted to the unit is brass C.P.

A wall- mounting clamp is provided with the unit and is supplied with the rubber tubing coil of 5 meters.

S.S.-304 grade Argon Arc Welded joints.

Supplied with cord and plug to work on 220/230 volts A.C. supply.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003