

Product Name :
Naugra Refreshing Cellar

Product Code :
NLAB-TECHNICALAB23002

Description :

Naugra Refreshing Cellar

Technical Specification :

Refreshing Cellar

Ventilated cooling +2°C / +10°C.

Frame and tank made of stainless steel 18/10 and cellar fully equipped with shelves.

Height of door: 1810 mm.

Soil isolated 10 mm.

Refrigeration group R404a.

Anti corrosion protection of the evaporator.

Defrosting by hot gas.

Re-evaporation of the condensates.

On / Off control switch, digital setting and display of the temperature and defrosting led.

Automatic defrosting and evaporation of water.

Ergonomic reversible door, with integral handfull and automatic back at 90°.

Insulation thickness 80 mm. 0% CFC.

Joint : 4 magnetic faces

Options :

Trolley with guides 13 floors for grids 600 x 400, réf. 050HAA.

Trolley with guides 13 floors for grids 630 x 510, réf. 050HBA.

Rolling platform for the transportation of the grids 630 x 510, réf. 050FCA.

Essential requirements:

Power of the group 450 W.

Maximal absorbed power 500 W.

Power supply : 230 V 50 Hz, single phase.

Interior volume designed for a trolley with guides (in option): (L x p x h) 705 x 740 x 1810 mm.

External dimensions: (L x W x H) 1045 x 1060 x 2300 mm.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003