

Product Name :
Nauga Semi Automatic Thermoscealing Machine

Product Code :
NLAB-TECHNICALAB23004



Description :

Nauga Semi Automatic Thermoscealing Machine

Technical Specification :

Semi Automatic Thermoscealing Machine

Semi Automatic Thermoscealing Machine cycle.

Down motion, pressing time and jack upward motion controlled by electronic timer.

Control of automatic cycle by positioning of pot support.

Precise adjustment of temperature done by probe and regulator.

Essential requirements

Air consumption : 0,5 m³ / hour- pressure 6 bars (compressed air not supply).

Single phase power : 220 Volts – 750 Watts.

Rhythm of working : 250 to 350 pots / hour.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003