


Product Name :
Nauga Single Disc Polishing Machine with controller

Product Code :
NLAB-ENGINEERINGLB84004



Description :

Nauga Single Disc Polishing Machine with controller

Technical Specification :

Metallography specimen polishing machine is suitable for the polishing of the specimen which has been grinded.

Single Disc Polishing - 8" Disc

The specimen surface is very smooth after processing and can be used to observe and measure the metallography structure of specimen under the microscope.

Table top design with continuously variable speed.

Quick change disc design enables discs with different emery grits/polish cloths to be used in desired steps.

Specification:

1/2HP high torque DC motor.

8" disc diameter Standard.

Thyristor controlled speed controller.

10" disc diameter option.

Many discs can be used.

Digital speed indicator - 50 to 1200 rpm.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003