

Product Name :
Naugra Sling Psychrometer Kit

Product Code :
NLAB-MLOGYLAB290006



Description :

Naugra Sling Psychrometer Kit

Technical Specification :

Measure relative humidity with this simple sling psychrometer. Consists of wet and dry bulbs and Celsius thermometers that rotate on a plastic handle. Correlate with other instruments in weather investigations.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003