

Product Name :
Nauga Specimen Mounting Machine

Product Code :
NLAB-ENGINEERINGLB30004



Description :

Nauga Specimen Mounting Machine

Technical Specification :

With firm commitment to quality, we are presenting our clients a vast range of **Specimen Mounting Machine**. These offered mounting machine is manufactured by the experts employing advanced technology to ensure its quality. Our esteemed clients can avail these mounting machines in various specifications according to their needs. Further, the offered range can be availed in different models and designs as per requirements of our clients at market leading prices.

Features:

- Simple construction
- Functional efficiency
- Resistance against corrosion

Stainless steel
Voltage
Digital
Electric Source
Power

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003