

Product Name :
Nauga Torque Measurement Trainer

Product Code :
NLAB-ENGINEERINGLB15002



Description :

Nauga Torque Measurement Trainer

Technical Specification :

Description:

Torque Lever

Length: - 0.5 Meter

Material: - Mild Steel

Digital Indicator: -

Least Count: - 0.1 Kg-M

Digit: - 3 ½

Weights: - Max 2 kg in steps of 100 gm

Powder coated enclosure.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003