

**Product Name :**  
Naugra Yoghurts Drier 140 Pots

**Product Code :**  
NLAB-TECHNICALAB23003

**Description :**

Naugra Yoghurts Drier 140 Pots

**Technical Specification :**

Yoghurts Drier 140 Pots

Capacity : 140 pots

Enamelled body drier with reinforced insulation

Adjustable and uniform temperature by fan

Can receive cardboard, plastic or glass pot.

Non reversible door (opening on the right when you look at the front view)

Can receive till 4 removable shelf wicker rack Dim. : 425 x 400 mm

**Dimensions and weight**

External dimensions : Height : 870 mm Width : 600 mm Depth : 520 mm

Internal dimensions : Height : 600 mm Width : 385 (430) mm Depth : 380 mm

Weight : 65 kg

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

