

Product Name :
Naugra Yoghurts Drier 140 Pots

Product Code :
NLAB-TECHNICALAB23003

Description :

Naugra Yoghurts Drier 140 Pots

Technical Specification :

Yoghurts Drier 140 Pots

Capacity : 140 pots

Enamelled body drier with reinforced insulation

Adjustable and uniform temperature by fan

Can receive cardboard, plastic or glass pot.

Non reversible door (opening on the right when you look at the front view)

Can receive till 4 removable shelf wicker rack Dim. : 425 x 400 mm

Dimensions and weight

External dimensions : Height : 870 mm Width : 600 mm Depth : 520 mm

Internal dimensions : Height : 600 mm Width : 385 (430) mm Depth : 380 mm

Weight : 65 kg

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

