

**Product Name :**  
Needle Penetrometer

**Product Code :**  
NLAB-EDUCATIONTVET11002



**Description :**

Needle Penetrometer

**Technical Specification :**

Complies with NF T66-004

Analog display on Ø150 mm dial graduated 360 °, 0.1 mm division

Micrometric descent of the needle

Supplied with needle, weight 50 and 100 g and brass cups Ø55x35 mm & Ø70x45 mm

Transfer container

Mirror for easy adjustment

IP 38C thermometer from +23 to + 26 ° C

Spare needle

Thermostatic bath

Closed circuit cooler

Sample cell with coil

Pipes and fittings

Hot air bath for softening samples before testing

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003