

**Product Name :**  
Nelaton Catheters

**Product Code :**  
NLABMEDICAL-M3030001



**Description :**

Nelaton Catheters

**Technical Specification :**

The Nelaton catheters are used for bladder catheterization for drainage purpose.

Soft, PVC material Nelaton catheter used for urine drainage.

Tube proximal end design with funnel connector to connect with urine bag.

Two lateral eyes are on tip for drainage without obstruction.

Radio opaque line provided through out the catheter for X-ray visualization.

Connector comes with different color for easy size identification of tube.

Sterile, single use.

Size: 6FG to 20FG, Length: 20 & 40 cms

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003