

Product Name :
Neonatal Intensive Care Incubator

Product Code :
NLABMEDICAL-M2740001



Description :

Neonatal Intensive Care Incubator

Technical Specification :

The Neonatal Intensive Care Incubator has been specially designed to provide the safest and most stable environment for the critical neonate.

With a wide range of accessories, it can be tailored to meet any situation arising in the care of the infant.

The excellent nursing facilities, safe and reliable operation, ease of cleaning and sturdy construction of the machine make it a vital part of the equipment in the hospital.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003