

Product Name : Neonatal Phototherapy Unit Cfl-A	Product Code : NLABMEDICAL-M2740003
Description :	
Neonatal Phototherapy Unit Cfl-A	
Technical Specification :	
<p>CFL module For the first time in the world, six lamps are used for speedy and effective treatment.</p> <p>The CFL's provide blue light at ~18 ?W/cm 2/nm. CFLs last five times longer than ordinary tube lights, requiring less frequent replacement.</p> <p>The CFL comes optionally with four blue and two white lamps for examination of the infant.</p> <p>Focused blue light</p> <p>Light from the CFL module is focused on the baby bed, avoiding scatter into the care-provider space.</p> <p>Heat dissipation</p> <p>Special vents are provided on the light unit to divert heat, avoiding undesirable warming of body.</p> <p>Compatibility with other equipment</p> <p>The tiltable light unit may be placed over any bed or incubator and in combination with a radiant warmer or another CFL unit for additional intensity.</p> <p>The height of the unit is adjustable, so that it is suitable for different bed heights.</p> <p>The T-Shaped base slides conveniently out of the way under other equipment.</p> <p>Accessory</p> <p>The CFL is provided with separately operated white lights.</p> <p>A chronometer is provided to monitor lamp usage.</p> <p>An additional chronometer may be used for management of treatment.</p>	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003