

Product Name : Neonatal Radiant Warmer-A	Product Code : NLABMEDICAL-M2740005
Description :	
Neonatal Radiant Warmer-A	
Technical Specification :	
IR radiation is beamed from a quartz heater by a parabolic reflector.	
This provides focused warmth rapidly on the desired area.	
The height of the IR unit may be adjusted for different bed heights.	
The low T-shaped base slides out of the way under other equipment.	
Infant temperature is measured quickly by a high sensitivity thermostat probe and displayed on a large, easy-to-read LED display.	
The control panel is ergonomically designed with soft-touch keys.	
The NWS is equipped with dazzle-free observation light for examination of the infant.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003