

Product Name :
Neonatal Radiant Warmer-A

Product Code :
NLABMEDICAL-M2740005

Description :

Neonatal Radiant Warmer-A

Technical Specification :

IR radiation is beamed from a quartz heater by a parabolic reflector.

This provides focused warmth rapidly on the desired area.

The height of the IR unit may be adjusted for different bed heights.

The low T-shaped base slides out of the way under other equipment.

Infant temperature is measured quickly by a high sensitivity thermostat probe and displayed on a large, easy-to-read LED display.

The control panel is ergonomically designed with soft-touch keys.

The NWS is equipped with dazzle-free observation light for examination of the infant.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003