

Product Name :
Neonate Scalp Venipuncture Model

Product Code :
NLAB-MEDICALMODELAB11040



Description :

Neonate Scalp Venipuncture Model

Technical Specification :

This model has main veins of scalp, including venae temporales superficals vein and venae frontales, and can be used to infusion, injection and phlebotomy.

There is obvious feeling when needle penetrates into the veins.

Blood flash back indicates proper insertion;

These operations won't cause leaks.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003