

Product Name : Neuro Attachment	Product Code : NLABMEDICAL-M2760005
---	---

--

Description :

Neuro Attachment

Technical Specification :

We offer a supreme range of Neuro Attachment, which is suitable for brain surgery.

This Neuro attachment has been developed & manufactured with innovative quality to meet the demanding needs of Neurosurgeons.

Key Structure

Provide 360-degree rotation positioning of the head system.

Adopts all major headrest and skull clamps to the Mayfield three head pin base unit.

Exclusive spring loaded mechanism provides assured stability and a secured lock.

The skull clamp has the additional feature of a swivelling c-arm.

This swiveling capacity greatly facilitates positioning the pins on the head around nerves, vessels. and thin bone in the temporal region.

For better result, the arm is locked by control knob and the width is adjusted by a ratchet mechanism.

The three cone-shaped pins available in adult and child sizes are interchangeable and may be removed for sterilization.

This attachment, originally designed for transsphenoidal hypophysectomy is used for all prone position procedures.

It is attached to the adjustable basic unit by means of transitional member and the gardner clamp adaptor.

It is padded with soft pliable sponge rubber for equalized weight distribution and incorporates the new quick-change pads.

The horse shoe width can be quickly adjusted to accommodate various head sizes.

A removable, adjustable pulley bar is included if traction is desired.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003